










Open MRI Comparison Chart

SIZE / POWER GUIDE

Oregon Advanced Imaging has created this MRI comparison chart to help medical offices understand the Open MRI technology available in the Rogue Valley. All specs have been provided by the manufactures.

			<p>Distance from face to scanner: 12"</p> <p>Image quality/strength: 1.5 tesla</p> <p>Max weight: 550 lbs.</p> <p>Opening height: 22"</p> <p>Opening width: 28</p> <p>Field of view: 19.7"</p> <p>Over 60% of scans can be completed with patients head out of scanner.</p>
<p>Oregon Advanced Imaging</p>		<p>Espree (1.5 tesla)</p>	<p>Specs</p>
			<p>Distance from face to scanner: 7.5"</p> <p>Image quality/strength: 1.2 tesla</p> <p>Max weight: 500 lbs.</p> <p>Opening height: 17.5"</p> <p>Opening width: 32.3"</p> <p>Field of view: 17.5"</p>
<p>Medford Open MRI</p>		<p>Oasis (1.2 tesla)</p>	<p>Specs</p>
			<p>Distance from face to scanner: 8.5"</p> <p>Image quality/strength: 1.0 tesla</p> <p>Max weight: 550 lbs.</p> <p>Opening height: 18.5"</p> <p>Opening width: 31.5"</p> <p>Field of view: 17.5"</p>
<p>Southern Oregon Open MRI</p>		<p>Panorama (1.0 tesla)</p>	<p>Specs</p>